

# BESA License Key 3.0 - Instructions for Use

Documents are subject to workflows, and content shall only be used if the content is approved (indicated by a green ring-symbol)

 Approved, Changed on: February 5th 2026 at 2:15:24 PM (UTC)



---

## Instructions for Use

---

### Purpose

This document contains the instructions for use and safety instructions, as well as technical description for this BESA License Key 3.0 in English language.

### Instructions for Use

The referenced term “BESA License Key 3.0” includes all sub-versions of the product.

This document, and following documents form the complete Instructions for Use for this BESA product:

- BESA License Key 3.0 - Intended Use and Product Classification
- BESA License Key 3.0 - Instructions for Use
  - How to update BESA Licenses
  - How to use BESA licenses on a local network
  - Communication problems with Sentinel License Managers
- BESA License Key 3.0 - External Release Notes
- BESA License Key 3.0 - Installation and Un-Installation
- BESA License Key 3.0 - Product Symbols Meaning

Current and previous versions of documents are available on BESA website:



<https://www.besa.de/downloads/user-documentation/>

It is mandatory to read the complete Instructions for Use documents, before using the product.

At any time during the lifetime of this product, you can request the above-mentioned documents in printed paper form at no additional cost.

## Product Name

BESA License Key 3.0

The referenced term “BESA License Key 3.0” includes all subversions of the product. BESA License Key 3.0 is an accessory product to BESA medical devices.

## Product Picture



## Product UDI

<b>UDI</b>	
UDI-DI:	(1)4262364230059
UDI-PI:	NOTE: the UDI-PI is individually for each License Key 3.0, including the serial number. It is available on the label applied to BESA Box.

## Manufacturer

	BESA GmbH Freihamer Strasse 18 82166 Gräfelfing – Germany
---	---

## Symbols used for Product

<b>CE</b>	A CE Declaration of Conformity under non-medical EU legislation, it covers: EMC Directive 2014/30/EU RoHS 2011/65/EU Electrical / ICT safety (EN IEC 62368-1) Issued by Thales DIS France SAS Due to the size of the BESA License Key, it is not marked on the device.
-----------	---

Further symbols used are available in document “BESA License Key 3.0 - Product Symbols Meaning”.

## Life cycle

There is no end of life for this version of the BESA License Key.

## Incident

In case of a serious incident, an incident or near-incident caused by the BESA License Key, the licensee should directly report to [incident@besa.de](mailto:incident@besa.de)<sup>1</sup>, follow the guidelines stipulated on the BESA website. (<https://www.besa.de/contact/incident/>).

## Safety Instructions and Technical Description

### Licensing System

BESA products are protected by the Sentinel LDK licensing system as developed by Thales.

The licensing system requires a USB key – that is, the BESA License Key 3.0 – and the device drivers to access the USB key.

If you selected the Install BESA License Key drivers option at the final stage of the installation, the device drivers for the BESA License Key 3.0 have been installed automatically on your computer. In order for BESA products to be run, a BESA License Key with valid license information must be attached to a USB port on a computer in your local area network. If you have not yet received it, you can obtain the BESA License Key from our sales department by contacting [sales@besa.de](mailto:sales@besa.de)<sup>2</sup>.

If you have received a license update for your BESA License Key and you are not sure how to apply it, please refer to “[How to update BESA licenses](#)”<sup>3</sup>.

---

1. <mailto:incident@besa.de>  
2. <mailto:sales@besa.de>  
3. <https://my.hidrive.com/lnk/4JnJV23>

If you have attached the BESA License Key 3.0 (USB key) directly to your computer or server and are not able to start the relevant BESA product, please refer to “[How to use BESA licenses on a local network](#)”<sup>4</sup> or “[Communication problems with Sentinel License Managers](#)”<sup>5</sup>.

## Serial Number / Dongle ID

If you want to know the serial number of your BESA License Key, please make sure that BESA License Key drivers (HASPUseSetup.exe, <https://www.besa.de/home/downloads/besa-license-key/>) are installed.

Connect the BESA License Key to your local PC and open the link <http://localhost:1947/><sup>6</sup> in any internet explorer of your choice.

The serial number can be found in the shown table, in the “Key ID” column (the key with the type “Sentinel HL Max” and Vendor ID “YWCXE (87128)”). If more than one BESA License Key is connected to your local PC, you can use the “Blink on” button: the corresponding key starts to blink.

## System Requirements

The BESA License Key 3.0 is an essential accessory product to BESA products, and shall either be subject to constant supervision, or shall be connected to machines located in areas with access restriction. Failure of the IT-network and/or operating system on which the BESA License Key 3.0 is used, or failure of the BESA License Key 3.0 in general, may render BESA products which use the BESA License Key 3.0 as an accessory as unusable.

- For the purposes of communication between BESA products and the BESA License Key 3.0, the manufacturer uses TCP and UDP ports 1947. Please make sure to open these ports within your network or IT environment in order to enable the required communication. Once the license key runtime environment is installed, you can find more information on settings, security and general topics concerning the BESA License Key 3.0 under the following link:

<http://localhost:1947> > Help

- HASPUseSetup.exe installs the runtime environment and drivers required by the BESA License Key 3.0. The installed version is listed under “Run-time Installer” under the following link:

<http://localhost:1947> > Diagnostics

Important: Make sure to regularly check website <https://www.besa.de/home/downloads/besa-license-key/> and apply updates for HASPUseSetup.exe if available. The supplied runtime environment is provided by BESA “as is”, and according to the manufacturer is compatible with computers running operating systems as listed on the BESA website.

---

4. <https://my.hidrive.com/lnk/unnlJlhi>

5. <https://my.hidrive.com/lnk/d0nop5IK>

6. [http://localhost:1947/\\_int\\_/devices.html](http://localhost:1947/_int_/devices.html)

- The BESA License Key 3.0 is designed to be attached to a USB 2.0 or USB 3.0 port. Please note that you may need to use an adapter to attach the BESA License Key 3.0 to your PC (the adapter is not included).

Please note, that changes of the hardware and software platforms (e.g., with patches/ updates of antivirus/firewall software, system libraries, firmware, and others) and/or changes to the IT network may result in the possibility of safety and/or security vulnerabilities, and/or reduced functionality of BESA License Key 3.0. It is recommended to use the system combinations as defined within this document.

Please note, that the use of BESA License Key 3.0 on an IT-network could result in previously unidentified risks to patients, users or third parties, and that the user is responsible to identify, analyze, evaluate, and control these risks. Any subsequent changes to the IT-network could introduce new risks and require additional analysis.

## **Data Protection**

BESA License Key 3.0 does not store any data.

## **Use in combination with other products**

The license key is allowed to be used as an accessory product for the BESA products as defined in document Intended Use.

BESA releases are constantly checked with other installations (e.g. MS Office etc.) on the same hardware or software platforms as defined in section System Requirements.

If any incompatibilities are observed, they are documented within the External Release Notes of the product. For further questions and information please use the channels provided on <http://www.besa.de>.

## **Application Environment**

Please note: It is the user's responsibility to ensure the safety of combined medical or medical-&-non-medical-devices particularly installed in the patient environment (whether in institutes, clinics or practitioners location).

The user of our product has to ensure that such combination fulfills the same safety condition as a single medical device according to IEC 60601-1.

Such combination with any non-medical device establishes an "in-house medical system" – observe the latest Medical Device Operator Ordinance (e.g. MedBetreiberVerordnung, Germany 2020) or your national Operator Laws. We recommend you to document / retrace our software version/updates installed in your facility / network in a MPG-Medizingeräte-Buch (medical facilities book), similar to those needed for hardware.

This BESA product does not affect the possibility of using any medical computer within patient environment.

BESA GmbH is not liable for any use/installation beyond the defined PC network.

## **Installation and Un-Installation**

Please refer to a separate document "[BESA License Key 3.0 - Installation and Un-Installation](#)<sup>7</sup>".

## **Exclusions for Use**

N/A

## **Contraindication to patients**

N/A

BESA is not liable for the use outside of the rules defined in the instructions for use documents.

## **Manufacture, Transport, Storage and Usage conditions**

### **Environmental conditions for Transport, Storage and Usage**

During extended periods of transport, storage and usage please observe the following conditions:

Temperature: -40°C to 85°C

Humidity: 0% to 95%

### **Disposal information**

The BESA License Key 3.0 and the box in which it is delivered must be disposed of in compliance with national guidelines on environmental protection.

### **Countries and Markets**

BESA License Key 3.0 usage is related to BESA products and the countries and markets defined for these products.

## **Page ID**

243007533

---

7. <https://my.hidrive.com/lnk/43noJydP>

## Change History (valid revision number)

Version	Effective Date	Approvals	Version Notes
1	Thu, Feb 22, 2024, 09:24:51 AM	>> *Action: approve, MK, Approval: QM Approval (State: QM Approval, Thu, Feb 22, 2024, 09:24:51 AM)	None
2	Wed, Apr 10, 2024, 02:17:39 PM	>> *Action: approve, MK, Approval: QM Approval (State: QM Approval, Wed, Apr 10, 2024, 02:17:38 PM)	Added information on how to get the serial number of the BESA License Key in ACC
3	Thu, Feb 5, 2026, 02:15:24 PM	>> *Action: approve, OK, Approval: Approved (State: Non-QM Approval, Thu, Feb 5, 2026, 02:15:24 PM)	None

## Attachments

Attachment	Description
N/A	N/A

## Template Change History

TEMPLATE	TEMPLATE Version	Effective Date	Version Notes
Instructions for Use Template	1	Tue, 14 Mar 2023 08:34:52 GMT	None
Instructions for Use Template	2	Tue, 21 Mar 2023 14:06:35 GMT	None
Instructions for Use Template	3	Tue, 02 May 2023 08:19:51 GMT	None
Instructions for Use Template	4	Thu, 25 May 2023 06:27:29 GMT	None

TEMPLATE	TEMPLATE Version	Effective Date	Version Notes
Instructions for Use Template	5	Mon, 26 Jun 2023 14:33:24 GMT	None

---

**Printed version is NO controlled document and NOT subject to modifications**