

BESA License Key - Update History

Sentinel HL NetTime 10 (V. 6.08):

Info

The Sentinel HL NetTime 10 keys serve as copy protection for BESA software and can protect one or multiple programs.

Individual features of an application can be activated with so-called HASP Features (slots) within the key.

These USB keys provide the license/s locally or within a network. The maximum number of concurrent users for each HASP feature is ten.

Improvements

Updated versions V. 4.31/V. 4.34 to 6.08:

BESA License Keys were updated with an alternate smart card chip. These keys are fully backward compatible. To help visually identify the different key version: The silver keyring was replaced with a gold keyring. (Please refer to roxtra, document ID 92388, Rev. 001/09.2021)

Bug fixes

None

Known Issues

None

HASPUserSetup.exe (V. 8.21):

Info

HASPUserSetup.exe is a GUI-based installation program to install the run-time environment for the BESA License Keys.

This utility is supplied by BESA license key manufacturer Thales and is provided by BESA as is.

The utility and the corresponding readme file are located in **Utilities\License\System** folder of the installation directory of BESA programs.

Improvements / Bug fixes / Known Issues

Please refer to the **Readme.html** file in **Utilities\License\System** folder of the installation directory of BESA programs.

BESA License Key - Update History



November 2021

Additional Known Issues (older versions)

<i>Affected Tool + Version</i>	<i>Publication Date</i>	<i>Information</i>
HASPUserSetup.exe (V. 7.100 and older)	November 2021	[KB0021742] Upgrade of Windows 10 to version 2004 may fail when working with run-time environment installed by HASPUserSetup.exe version 7.100 or earlier Please refer to the original manufacturer article for more information and workarounds: link to KB0021742