



Product Symbols



Product Symbols



The CE marking certifies that this product fulfills the essential requirements of the Medical Devices Directive MDD 93/42/EEC.

The number 0197 represents the identification number of the Notified Body.

Symbol significant for:

BESA Epilepsy, BESA MRI, BESA Research



The CE marking certifies that this product fulfills the essential requirements of the Medical Devices Directive MDD 93/42/EEC.

Symbol significant for:

**BESA Connectivity, BESA Statistics,
BESA License Key**

Product Symbols



Patient environment – critical area for interactions of patient, applied parts, electronic medical devices – special medical PC / devices to use if the product is installed on a PC or notebook within the patient environment.



Symbol significant for:

**BESA Connectivity, BESA Epilepsy,
BESA MRI, BESA Research, BESA Statistics,
BESA License Key**



Identifies the manufacturer of the product.

Symbol significant for:

**BESA Connectivity, BESA Epilepsy, BESA MRI,
BESA Research, BESA Statistics**



Indicates that it is mandatory to read and to follow instructions for use before using product.

Symbol significant for:

**BESA Epilepsy, BESA MRI, BESA Research,
BESA Statistics, BESA License Key**

Product Symbols



Shows temperature limits for storage, transport and usage.

Symbol significant for:

BESA Connectivity, BESA Epilepsy, BESA MRI, BESA Research, BESA Statistics, BESA License Key



Shows humidity limits for storage, transport and usage.

Symbol significant for:

BESA Connectivity, BESA Epilepsy, BESA MRI, BESA Research, BESA Statistics, BESA License Key



FCC-marking - Compliance to US requirements with respect to electromagnetic interference (EMC) from the device.

Symbol significant for:

BESA License Key

Product Symbols



The WEEE directive sets collection, recycling and recovery targets for all types of electrical goods with the purpose of dramatically reducing the amount of harmful waste materials that are placed into landfills each year.

Symbol significant for:
BESA License Key



Chinese eMark (environmental Marking - covering the restricted use of hazardous substances in electronic and electrical devices). This product contains at least one hazardous substance but can be used safely during its environmental protection use period (EPUP), which is as indicated by the number in the center. After its defined EPUP the product should be recycled or disposed of properly.

Symbol significant for:
BESA License Key

Product Symbols



Indicates the date after which the medical device is not to be used.

For BESA products only related to Chinese eMark on BESA License Key.

Symbol significant for:
BESA License Key



Indicates the Serial Number of the device.

Symbol significant for:
BESA License Key



Indicates that the device is a Medical Device.

Symbol significant for:
**BESA Connectivity, BESA Epilepsy, BESA MRI,
BESA Research, BESA Statistics**

Product Symbols



Indicates the Unique Device Identifier of the device.

Symbol significant for:

**BESA Connectivity, BESA Epilepsy,
BESA MRI, BESA Research, BESA Statistics,
BESA License Key**



Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.

Symbol significant for:

BESA MRI

BESA GmbH
Freihamer Str. 18
82166 Graefelfing – Germany

Phone +49. 89. 89809966
E-mail info@besa.de
Web www.besa.de

© Copyright 2020

Document name BESA_Booklet_Product_Symbols
Revision number 004/06.2020
Document ID 86227

