





















| | | | | |
|---------------------------------------|--------------------------------|---|--|--|
| BESA Product Symbols - Meaning | |  BESA® | | |
| Document ID 93543 | Revision Number 004/07.2022 | Edited by Kornwebel, Michael 2022-07-05 | Reviewed by: Scherg, Tobias 2022-07-05 | Released by / Date Scherg, Tobias 2022-07-05 |


| Symbol | Meaning | Applicable for BESA Products |
|---|---|--|
|  | The CE marking certifies that this product fulfills the essential requirements of the Medical Devices Directive MDD 93 /42 / EEC. The number 0197 represents the identification number of the Notified Body | BESA Research 7.1.2.1 BESA MRI 3.0.0.0 |
|  | The CE marking certifies that this product fulfills the essential requirements of: the Medical Devices Directive MDD 93 /42 / EEC; or the Medical Device Regulation MDR (EU) 2017/745 | BESA Connectivity 1.0.1.0 BESA License Key 2.0.0.0 |
|  | Indicates the medical device manufacturer | BESA Research 7.1.2.1 BESA MRI 3.0.0.0 BESA Connectivity 1.0.1.0 BESA License Key 2.0.0.0 |
|  | Indicates the need for the user to consult the instructions for use | BESA Research 7.1.2.1 BESA MRI 3.0.0.0 BESA Connectivity 1.0.1.0 BESA License Key 2.0.0.0 |
|  | Indicates the item is a medical device | BESA Research 7.1.2.1 BESA MRI 3.0.0.0 BESA Connectivity 1.0.1.0 |
|  | Indicates the manufacturer's serial number so that a specific medical device can be identified | BESA License Key 2.0.0.0 |
|  | Indicates a carrier that contains unique device identifier information | BESA Research 7.1.2.1 BESA MRI 3.0.0.0 BESA Connectivity 1.0.1.0 BESA License Key 2.0.0.0 |
|  | Indicates that the original medical device information has undergone a translation which supplements or replaces the original information | BESA Research 7.1.2.1 BESA MRI 3.0.0.0 BESA Connectivity 1.0.1.0 BESA License Key 2.0.0.0 |
|  | Indicates the temperature limits to which the medical device can be safely exposed | BESA License Key 2.0.0.0 |

| | |
|--|---------------------------------------|
| Template document name: BESA Product Symbols-Meaning | Template edited on/by: 2022-Q2/ MK |
| Template valid from: 2022-05-09 | Template released on/by: 2022-05 / MK |
| Template Document ID: 93543 | |

| | | | | |
|---------------------------------------|--------------------------------|---|--|--|
| BESA Product Symbols - Meaning | |  BESA® | | |
| Document ID 93543 | Revision Number 004/07.2022 | Edited by Kornweibel, Michael 2022-07-05 | Reviewed by: Scherg, Tobias 2022-07-05 | Released by / Date Scherg, Tobias 2022-07-05 |

| Symbol | Meaning | Applicable for BESA Products |
|---|---|------------------------------|
|  | Indicates the range of humidity to which the medical device can be safely exposed | BESA License Key 2.0.0.0 |
|  | Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information | BESA License Key 2.0.0.0 |
|  | FCC-marking - Compliance to US requirements with respect to electromagnetic interference (EMC) from the device | BESA License Key 2.0.0.0 |
|  | This symbol indicates that when the end-user wishes to discard this product, it must be sent to separate collection facilities for recovery and recycling. By separating this product from other household-type waste, the volume of waste sent to incinerators or land-fills will be reduced and natural resources will thus be conserved. | BESA License Key 2.0.0.0 |
|  | Defines the Environment-Friendly Use Period (EFUP) as 10 years after production. Note: production date of respective BESA License Key is included in UDI info, in format YYMMDD, following application identifier: (11). | BESA License Key 2.0.0.0 |
|  | Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences | BESA MRI 3.0.0.0 |
|  | Refer to instruction manual/booklet | BESA License Key 2.0.0.0 |
|  | Product certification which ensures the conformity of products to Korean safety standards | BESA License Key 2.0.0.0 |
|  | Symbol refers to BESA Research | BESA Research 7.1.2.1 |
|  | Symbol refers to BESA MRI | BESA MRI 3.0.0.0 |


| | |
|--|---------------------------------------|
| Template document name: BESA Product Symbols-Meaning | Template edited on/by: 2022-Q2/ MK |
| Template valid from: 2022-05-09 | Template released on/by: 2022-05 / MK |
| Template Document ID: 93543 | |

| | | | | |
|---------------------------------------|-----------------|---|------------------------------|------------------------------|
| BESA Product Symbols - Meaning | |  BESA [®] | | |
| Document ID | Revision Number | Edited by | Reviewed by: | Released by / Date |
| 93543 | 004/07.2022 | Kornweibel, Michael 2022-07-05 | Scherg, Tobias 2022-07-05 | Scherg, Tobias 2022-07-05 |

| Symbol | Meaning | Applicable for BESA Products |
|---|---|---|
|  | Symbol referring to, and used with BESA Connectivity | BESA Connectivity 1.0.1.0 |
|  | Symbol refers to BESA Statistics | BESA Statistics 2.1.0.0 |
|  | Symbol refers to BESA pipeline | Not a BESA product. |
|  | Symbol used with BESA MRI | BESA MRI 3.0.0.0 |
|  | Symbol used with BESA Research | BESA Research 7.1.2.1 |
|  | Symbol used with BESA Statistics | BESA Statistics 2.1.0.0 |
|  | Symbol used with BESA License Search Note: support application to BESA programs | BESA Research 7.1.2.1 BESA MRI 3.0.0.0 BESA Connectivity 1.0.1.0 BESA Statistics 2.1.0.0 |
|  | Symbol used with ConfigureBesaMatlabInterface, SetSecAttrToFullCtrlBesaResearch, UpdateIni InstallReader (reader: Multisphere) Note: used with the BESA Research modules | BESA Research 7.1.2.1 |
|  | Symbol used with Diagnostics Note: support application to BESA programs | BESA Research 7.1.2.1 BESA MRI 3.0.0.0 BESA Connectivity 1.0.1.0 |
|  | Symbol used with InstallReader (different reader) Note: used with the BESA Research module | BESA Research 7.1.2.1 |
|  | Symbol used with InstallReader (different reader) Note: used with the BESA Research module | BESA Research 7.1.2.1 |

| | |
|--|---------------------------------------|
| Template document name: BESA Product Symbols-Meaning | Template edited on/by: 2022-Q2/ MK |
| Template valid from: 2022-05-09 | Template released on/by: 2022-05 / MK |
| Template Document ID: 93543 | |

| | | | | |
|---------------------------------------|--------------------------------|--|--|--|
| BESA Product Symbols - Meaning | |  | | |
| Document ID 93543 | Revision Number 004/07.2022 | Edited by Kornwebel, Michael 2022-07-05 | Reviewed by: Scherg, Tobias 2022-07-05 | Released by / Date Scherg, Tobias 2022-07-05 |

| Symbol | Meaning | Applicable for BESA Products |
|---|---|--|
|  | Symbol used with ACC Admin Control Center | BESA Research 7.1.2.1 BESA MRI 3.0.0.0 BESA Connectivity 1.0.1.0 BESA License Key 2.0.0.0 |

| | |
|--|---------------------------------------|
| Template document name: BESA Product Symbols-Meaning | Template edited on/by: 2022-Q2/ MK |
| Template valid from: 2022-05-09 | Template released on/by: 2022-05 / MK |
| Template Document ID: 93543 | |