


BESA Research 7.1 - Product Symbols Meaning

Documents are subject to workflows, and content shall only be used if the content is approved (indicated by a green ring-symbol)

 Approved , Changed on: February 9th 2024 at 1:30:34 PM (UTC)



Symbols-Meaning

This document provides all symbols used with respect to:

BESA Research 7.1.3.0

UDI (human readable) related to this product is available in [BESA Research 7.1 - Instructions for Use](#)


AIDC (Barcode) UDI Labeling of Product Version:

BESA Research 7.1 is a software only product. It will be available for users by download from BESA website, no packaging, shipping involved. According to MDR Annex I, Chapter III section 23.1 (c), Article we will provide the human readable version of the UDI in the About box only.

Symbol	Occurrence at BESA Product or accompanying documents, website	Meaning	Source of Symbol
---------------	--	----------------	-------------------------

About box:

About BESA Research ✕



BESA® Research

MEG + EEG 7.1.3.0, January 2024

Copyright © 2024, BESA GmbH. All rights reserved!

The BESA Research software is licensed for use on a single computer system and is provided "as-is". All warranties of any kind are hereby disclaimed, including all warranties of fitness for a particular purpose.

In countries where the CE certification is not recognized, it is not allowed to use this software directly or indirectly for medical diagnosis and/or treatment of humans.
Not for diagnostics use in the United States of America!

Please consult the Getting Started guide and the Intended Use document before first use. A comprehensive help documentation is available from the Help menu.

BESA Research makes use of unmodified third-party libraries. Please see \Utilities\License\BESA Research - License Acknowledgements.pdf for license acknowledgements.


If you accept the above licensing condition, click 'Accept', else click 'Cancel'


CE0197



MD

UDI


(1)4262384230004
(8012)07010300

 BESA GmbH, Freihamer Str. 18,
82166 Graefelfing, Germany

 <https://www.besa.de/home/downloads/besa-research/>

CE0197	About Box	The CE marking certifies that this product fulfills the essential requirements of the Medical Devices Directive MDD 93 /42 / EEC. The number 0197 represents the identification number of the Notified Body	Use of CE sign regulated by Medical Devices Directive MDD 93 /42 / EEC; or the Medical Device Regulation MDR (EU) 2017/745
	About Box	Indicates the medical device manufacturer	ISO 15223-1:2021
 eIFU indicator	About Box	Indicates the need for the user to consult the instructions for use	ISO 15223-1:2021
<div style="border: 1px solid black; padding: 5px; display: inline-block;">MD</div>	About Box	Indicates the item is a medical device	ISO 15223-1:2021
<div style="border: 1px solid black; padding: 5px; display: inline-block;">UDI</div>	About Box	Indicates a carrier that contains unique device identifier information	ISO 15223-1:2021

	<p>Instructions for Use</p>	<p>Indicates that the original medical device information has undergone a translation which supplements or replaces the original information</p>	<p>ISO 15223-1:2021</p>
	<p>BESA Research, within software</p>	<p>Used as icon for: desktop shortcuts, Start Menu, Program info entry on target PC</p>	<p>372dpi</p>
	<p>BESA Box BESA Booklet Image Brochure</p>	<p>Symbol used in BESA booklet, refers to BESA Research</p>	<p>372dpi</p>
	<p>BESA Box BESA Booklet Image Brochure</p>	<p>Symbol used in BESA booklet, refers to BESA MRI</p>	<p>372dpi</p>
	<p>BESA Connectivity, within software BESA Box BESA Booklet Image Brochure</p>	<p>Used as icon for: desktop shortcuts, Start Menu, Program info entry on target PC Symbol used in BESA booklet, refers to BESA Connectivity</p>	<p>372dpi</p>
	<p>BESA Box BESA Booklet Image Brochure</p>	<p>Symbol used in BESA booklet, refers to BESA Statistics</p>	<p>372dpi</p>
	<p>BESA Box BESA Booklet Image Brochure</p>	<p>Symbol used in BESA booklet, refers to BESA pipeline</p>	<p>372dpi</p>
	<p>BESA MRI, within software</p>	<p>Used as icon for: desktop shortcuts, Start Menu, Program info entry on target PC</p>	<p>372dpi</p>
	<p>BESA License Search, within BESA software</p>	<p>Used by BESA programs. Used as icon for: Start Menu entry on target PC</p>	<p>372dpi</p>

	Tools, e.g. Configure BESA MATLAB interface, within BESA software	Used by BESA programs. Used as icon for: entry on target PC	372dpi
---	---	--	--------

Page ID

212369433

Change History (valid revision number)

Version	Effective Date	Approvals	Version Notes
1	Fri, Feb 9, 2024, 13:30:34 PM	>> *Action: approve, VR, Approval: Non-QM Approval (State: Non-QM Approval, Fri, Feb 9, 2024, 13:30:34 PM)	None

Template Change History

TEMPLATE	TEMPLATE Version	Effective Date	Version Notes
BESA Product Symbols Meaning	1	Thu, 16 Mar 2023 07:02:54 GMT	None
BESA Product Symbols Meaning	2	Tue, 21 Mar 2023 14:36:24 GMT	None
BESA Product Symbols Meaning	3	Tue, 02 May 2023 07:37:57 GMT	None
BESA Product Symbols Meaning	4	Tue, 02 May 2023 07:40:08 GMT	None
BESA Product Symbols Meaning	5	Wed, 24 May 2023 14:26:05 GMT	None
BESA Product Symbols Meaning	6	Fri, 23 Jun 2023 06:19:40 GMT	None
BESA Product Symbols Meaning	7	Mon, 09 Oct 2023 12:01:45 GMT	None

Printed version is NO controlled document and NOT subject to modifications