

BESA License Key 3.0 - Product Symbols Meaning

Documents are subject to workflows, and content shall only be used if the content is approved (indicated by a green ring-symbol)

 Approved, Changed on: February 17th 2026 at 9:20:13 AM (UTC)



Symbols-Meaning


This document provides all symbols used with respect to:

BESA License Key 3.0 (all sub-versions of the product included)

UDI (human readable) related to this product is available in “BESA License Key 3.0 - Instructions for Use”

AIDC (Barcode) UDI Labeling of Product Version:

AIDC (automatic identification and data capture) are created individually for each BESA License Key, and is available on the label of BESA Box. Example:



Symbol	Occurrence at BESA Product or accompanying documents, website	Meaning	Source of Symbol
	Available in accompanying document “BESA License Key 3.0 - Instructions for Use”	This product bears a CE sign related to non-medical EU legislation	CE mark as defined by the EU

	<p>Available in accompanying document "BESA License Key 3.0 - Instructions for Use"</p>	<p>Indicates the contact address for the manufacturer of the device</p>	<p>ISO 15223-1:2021</p>
	<p>Available on the label of packaging, BESA Box</p>	<p>Indicates the need for the user to consult the instructions for use</p>	<p>ISO 15223-1:2021</p>
	<p>Available on the label of packaging, BESA Box, and in accompanying document "BESA License Key 3.0 - Instructions for Use"</p>	<p>Indicates a carrier that contains unique device identifier information, is created based on the manufacturer's serial number so that a specific medical device can be identified</p>	<p>ISO 15223-1:2021</p>
	<p>BESA shipping packaging</p>	<p>Indicates the temperature limits to which the medical device can be safely exposed</p>	<p>ISO 15223-1:2021</p>
	<p>BESA shipping packaging</p>	<p>Indicates the range of humidity to which the medical device can be safely exposed</p>	<p>ISO 15223-1:2021</p>
	<p>BESA shipping packaging</p>	<p>Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information</p>	<p>ISO 15223-1:2021</p>

	<p>BESA shipping packaging</p>	<p>This symbol indicates that when the end-user wishes to discard this product, it must be sent to separate collection facilities for recovery and recycling. By separating this product from other household-type waste, the volume of waste sent to incinerators or landfills will be reduced and natural resources will thus be conserved.</p>	<p>Directive 2012/19/EU WEEE (Waste Electrical and Electronic Equipment)</p>
	<p>BESA License Key, marked on device</p>	<p>FCC-marking - Compliance to US requirements with respect to electromagnetic interference (EMC) from the device</p>	<p>Regulated by US, more info at: https://www.fcc.gov/ Note: provided to BESA License Key by Thales</p>
	<p>BESA License Key, marked on device</p>	<p>Product certification which ensures the conformity of products to Korean safety standards</p>	<p>Korean regulations</p>

Note: the declaration of conformity by manufacturer can be requested from BESA.

Page ID

242679869

Change History (valid revision number)

Version	Effective Date	Approvals	Version Notes
1	Mon, Mar 25, 2024, 03:23:33 PM	>> *Action: approve, OK, Approval: Approved (State: Non-QM Approval, Mon, Mar 25, 2024, 03:23:33 PM)	None
2	Tue, Feb 17, 2026, 09:20:13 AM	>> *Action: approve, MK, Approval: QM Approval (State: QM Approval, Tue, Feb 17, 2026, 09:20:13 AM)	None

Template Change History

TEMPLATE	TEMPLATE Version	Effective Date	Version Notes
BESA Product Symbols Meaning	1	Thu, 16 Mar 2023 07:02:54 GMT	None
BESA Product Symbols Meaning	2	Tue, 21 Mar 2023 14:36:24 GMT	None
BESA Product Symbols Meaning	3	Tue, 02 May 2023 07:37:57 GMT	None
BESA Product Symbols Meaning	4	Tue, 02 May 2023 07:40:08 GMT	None
BESA Product Symbols Meaning	5	Wed, 24 May 2023 14:26:05 GMT	None
BESA Product Symbols Meaning	6	Fri, 23 Jun 2023 06:19:40 GMT	None
BESA Product Symbols Meaning	7	Mon, 09 Oct 2023 12:01:45 GMT	None

Printed version is NO controlled document and NOT subject to modifications